ASSIGNMENT-1 (10 Marks)

- 1. Download a PDB file of your choice and read the remarks sections of the file. Mention the PDB id of the protein and describe your protein in a maximum of 2-3 lines
- 2. Is your PDB structure determined by NMR or X-RAY? Which technique can result in resolving hydrogen atoms and provide hydrogen co-ordinates?
- 3. Represent the protein in Licorice and color the protein according to its elements.
- 4. Does your protein have crystal waters? If yes, represent using VDW.
- 5. Change the background color to white, hide the axes and render the image in orthographic mode. (Use this layout for answering rest of the questions)
- 6. Show all hetero atoms (need not include waters) lying within 2.0 Angstrom of your protein in CPK form. List them.
- 7. Represent in VDW the first atom of your PDB ID using keyword: index and serial. What do your infer from the observation?
- 8. Represent the first 5 residues in CPK. Find the angle between any three non-consecutive atoms and display the angle value in red and atom labels in blue.
- 9. Load the same PDB twice and place them next to each other in different color.
- 10. Display the contact map of protein. What do you interpret from the graph?